

# NEXT RELEASES

## Release 2.08.00rc4

Pre-release available from [Abacus](#) (...).

Fixing two bugs introduced with 2.08.00rc3.

- 🐛1980 new ⇒ Full-text search in patient history not working
- 🐛1981 new ⇒ Auto-delete of uploaded files fails the second time

## Release 2.08.00rc3

Pre-release available from [Abacus](#) (15/12/2017).

- 🐛1890 tested OK ⇒ Improvements of Audit-Trail (Patient ID, incident Entries, ...)
- 🐛1945 tested OK ⇒ Attached files corrupted during import
- 🐛1959 tested OK ⇒ Gecamed Backup consumes DB connections until pool is depleted
  - 🐛817 closed ⇒ No ManagedConnections available within configured....
  - 🐛818 closed ⇒ Failing backups, JBoss refusing connection
- 🐛1961 fixed ⇒ File Lock preventing Movement of Temp Files during Patient Import
- 🐛1962 tested OK ⇒ Delete files after attaching
- 🐛1963 tested OK ⇒ Attached files cannot be modified, even though it looks that way
- 🐛1965 tested OK ⇒ Allow for creating incidents on past date while attaching files
- 🐛1968 tested OK ⇒ Replace underscore in patient IDs sent by the DICOM worklist server
- 🐛1972 tested OK ⇒ Automatic window adjustment broken in DICOM Viewer
- 🐛1975 tested OK (experimental feature) ⇒ Add Category to Incidents and Incident Entries - GUI

Seems to work in general (no obvious bugs, errors or broken features found). See my general remarks for the usage of this new feature in the ticket notes.

- 🐛1976 tested OK ⇒ Patient's invoice list not updated after automatic reload

One minor issue remains: invoice state icon not updated. Other fields are.

- 🐛1977 tested OK ⇒ Splash messages in GECAMed are shown too short.

## Release 2.08.00rc2

Pre-release available from our support partner [Abacus](#) (24/8/2017).

- 🐛1935 tested OK ⇒ Error when reloading the history in Form view
- 🐛1936 tested OK ⇒ Date default values in FormEditor don't work as expected
- 🐛1943 wontfix ⇒ Possible to print label with no patient loaded
- 🐛1944 tested OK ⇒ Default label template always remains visible even if custom templates defined

- 🐾 1945 tested OK ⇒ Attached files corrupted during import
- 🐾 1946 tested OK ⇒ Incomplete import of patient memos fixed
- 🐾 1947 tested OK ⇒ Incomplete import of invoice memo
- 🐾 1948 tested OK ⇒ Enhance debugging/logging output of the patient import
- 🐾 1950 tested OK ⇒ Add invoice number and date to “liste des actes facturés au patient”
- 🐾 1951 tested OK ⇒ Add Category to Incidents and Incident Entries.
- 🐾 1217 tested OK ⇒ Empty measurement entry created when importing patient from XML

## Release 2.08.00rc1

Pre-release available from our support partner [Abacus](#) (28/6/2017).

- 🐾 1396 fixed ⇒ Can't print invoice on freshly installed Linux Ubuntu system
- 🐾 1880 tested OK ⇒ Modify Invoice Templates to use Codes instead of hard linking Rates
- 🐾 1881 tested OK ⇒ Option to “serialize” multiple print copies
- 🐾 1883 tested OK ⇒ Enhance debugging/logging output of the ruleengine
- 🐾 1885 tested OK ⇒ Include global and user settings in PDF report (Help→ About→ Save info)
- 🐾 1888 tested OK ⇒ Arzt-Spalte bei Rechnungserfassung IMMER anzeigen
- 🐾 1889 tested OK ⇒ Remove / reduce JMS messages when changing width of waiting room
- 🐾 1890 tested OK ⇒ Improvements of Audit-Trail (Patient ID, incident Entries, ...)
- 🐾 1892 tested OK ⇒ Avoid creating Antecedents without shortcut
- 🐾 1896 tested OK ⇒ Add patients “storage” field to patient search
- 🐾 1897 tested OK ⇒ Appointment target date is not respected even though physician is available
- 🐾 1898 to be tested ⇒ Import maiden name from HL7 message
- 🐾 1901 tested OK ⇒ Warning when creating an invoice for non-treating physician
- 🐾 1902 tested OK ⇒ Remove partially printed bank account number in 2nd reminder letter
- 🐾 1905 reopened ⇒ Doppelte Zahlungen in GECAMed
- 🐾 1906 tested OK ⇒ Display Chapter of Rate in the Rate Edit Dialog.
- 🐾 1907 tested OK ⇒ Extend length of amount field for direct amount entry in invoices
- 🐾 1908 tested OK ⇒ Briefvorlage Kinesitherapie falsch (32 anstatt 64 séances max)
- 🐾 1910 tested OK ⇒ Display empty Incidents in the history and give the option to delete them.
- 🐾 1911 tested OK ⇒ GECAMed startup slow if Autocompleter has many entries.
- 🐾 1912 tested OK ⇒ Add Option to Print a Label on a Labelprinter
- 🐾 1719 tested OK ⇒ Print patient address on label printer
- 🐾 1913 tested OK ⇒ Error saving patient when changing health insurance
- 🐾 1914 to be tested ⇒ Client startup/login speed - Host identification
- 🐾 1915 tested OK ⇒ Improved prescription follow up by marking prescriptions that will receive results

See also in the User Manual [Attendre les résultats](#) and [Auf Ergebnisse warten](#)

- 🐾 1916 tested OK ⇒ Some values not readable in info\_document.pdf
- 🐾 1917 tested OK ⇒ Deleting isolated accidents results in invalid history entry results
- 🐾 1918 tested OK ⇒ Doctors email address on invoices
- 🐾 1919 tested OK ⇒ Bad dp patch version shown in the client log
- 🐾 1921 tested OK ⇒ Antecedents shortcut not stored if longer than 250 characters
- 🐾 1922 tested OK ⇒ Incomplete logging of appointment changes in agenda
- 🐾 1923 tested OK ⇒ Error when saving patient after prescription modification
- 🐾 1924 tested OK ⇒ Invoice status not updated if paid in the statement list

- 🐾 [1925](#) experimental feature ⇒ Implement DICOM Workflow and DICOM Modality Worklist Server
- 🐾 [1926](#) experimental feature ⇒ MWL query result incomplete if connected user is not a physician
- 🐾 [1928](#) tested OK ⇒ Saving an invoice on one client unsaves it on the other
- 🐾 [1929](#) won't fix ⇒ Can't change the amount of a tariff due to lack of privileges
- 🐾 [1930](#) tested OK ⇒ RCAM majoration not applied if invoice created by shortcut
- 🐾 [1931](#) tested OK ⇒ JMS exception thrown when closing / reopening a statement
- 🐾 [1934](#) fixed ⇒ Deletion of form entry not reflected in history
- 🐾 [1935](#) fixed ⇒ Error when reloading the history in Form view

From:

<https://gm.apps.lu/> - **GECAMed - Gestion de Cabinets Médicaux**

Permanent link:

[https://gm.apps.lu/releasenotes\\_nextrelease?rev=1526021206](https://gm.apps.lu/releasenotes_nextrelease?rev=1526021206)



Last update: **2019/12/09 10:20**